dPROG

Ghid de Testare

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
| Phase | Proestrus | Pre- LH surge | LH surge | Preovulation | At or near Ovulation | Diestrus |
|-----------------------|--------------------------------------|---|------------------------|---------------------------------------|--|---|
| PROG | | | | | | |
| ng/mL | <1.0 | 1.0 - 1.9 | 2.0 - 3.0 | 3.1 - 4.9 | 5.0 - 12.0 | |
| nmol/L | <3.2 | 3.2 - 6.3 | 6.4 - 9.6 | 9.6 - 15.8 | 15.9 - 38.2 | |
| Clinical Signs | Vaginal cytology cornification index | 80% - 100%* | 80% - 100%* | 80% - 100%* | 80% - 100%* | 0% - 80% |
| | Vaginal cytology | Serosanguinous discharge or erythrocytes in vaginal smear, Leukocytes None or few | | Leukocytes None or few** | Leukocytes None or few** | Leukocytes Few to many |
| | Vulvar turgor | Progressively increasing | Peak or just decreased | Start to decrease | Obviously decreased | Vulva soft but still enlarged |
| | Vaginal Mucosal Folds | White, slight shrinking and wrinkling | | White, grossly wrinkled and angulated | | White and pink or red, blotchy, nearly flat |
| Action | Retest in 3 - 4 days | Retest in 2 - 3 days | Retest in 1 - 2 days | Retest in 1 - 2 days | Retest in 1 day; determine breeding schedule.*** | Infertile or Rarely fertile |

* Typically 95 - 100% cornification by 2 - 8 days before LH surge.

** Typically none, as leukocytes are expected to be absent in smear from 4 - 10 days before the surge until 7 - 10 days after LH surge.

***Breeding schedule based on desired method of insemination and semen type.

The information contained herein is intended to provide general guidance only, individual bitches may vary. Any suggestion should not be the clinical judgement directly, please use in conjunction with other monitoring techniques, such as observation of clinical signs, vaginal cytology, ultrasound and additional hormonal testing.